


## Spare Parts List Incubator 8000 SC, US version




**Summary of pictures:**

Picture	Description
2 - 3	BASIC UNIT INCUBATOR 8000 SC
4	CARRIAGE WITH LIFTING COLUMN
5	PEDAL FOR HEIGHT ADJUSTMENT
6	PLATE FOR HEIGHT ADJUSTMENT
11	FAN WHEEL
12 - 15	HOUSING NC/SC/IC
16	HINGE BRACKET FOR HOOD
17	TAPPET FOR STANDARD HOOD
18	CATCH FOR SCAN. HOOD
19	SWIVEL CUPBOARD
20	MOUNT FOR CUPBOARD
23	ELECTRIC MAINS NC/SC, VARIABLE HEIGHT / FIXED HEIGHT
24	ELECTRIC MAINS NC/SC, VARIABLE HEIGHT
25	AGGREGATE
26	PCB SHEET METAL (PART OF AGGREGATE)
27	MOTOR
28	HUMIDIFIER
29	HANDLE
30	INCLINATION, RIGHT, LEFT
31 - 32	MODI. SET H2O/INC8000 SC TO IC
33	WATER CONTAINER SET
34	STANDARD HOOD W/O DOUBLEWALLS (OLD VERSION)
37	STANDARD HOOD W/O DOUBLEWALLS (OLD VERSION)
42	STANDARD HOOD DOUBLEWALLS (NEW VERSION)
35	FLAP FOR STANDARD HOOD
36	SET OF SWIVELLING WINDOWS
38	SCAND. HOOD W/O DOUBLEWALLS (OLD VERSION)
46	SCAND. HOOD W/O DOUBLEWALLS (OLD VERSION)
43	SCAND. HOOD DOUBLEWALLS (NEW VERSION)

Picture	Description
36	SET OF SWIVELLING WINDOWS
39	FLAP, BIG, FOR SCAND. HOOD
40 - 41	FLAP, SMALL, FOR SCAND. HOOD
47	TEMPERATURE-SENSOR
48	BORDERING.
49 - 55	PLUG-IN UNIT 8000 NC/SC
56	PCB DISPLAY AIR
57	PCB POWER PACK
58	PCB CPU INCUBATOR
59	PCB ANALOG, THERMON
60	PCB DISPLAY HUMIDITY
61	PCB DISPLAY SKIN
62	INSTALLATION SET RS 232 (OPTION)
63	CABLE HARNESS SKIN TEMPERATURE MEASUREMENT
66	SKIN TEMPERATURE SENSORS
67	MODIFICATION SET LANGUAGE
68	VENTILATION HOSE HOLDER
69	MODIFICATION SET USA
70	MODIFICATION KIT INCU. 8000 SC, USA
72	MODIFICATION KIT CORNER PIECE


## BASIC UNIT INCUBATOR 8000 SC

Picture 2



**BASIC UNIT INCUBATOR 8000 SC**

**Picture 3**




Item No.	Description	Part No.	Order-Code
1	MAINS CABLE UL		2M20453
3	FUSE LINK DIN 41662 T6,3A		1815172
4	AGGREGATE, <a href="#">see picture 25</a>		2M20615
5	SCREW AM 6X10 DIN 85-A2	1330764	
6-22	MODIFICATION KIT HOUSING SC/IC, <a href="#">see picture 12 - 15</a>	2M21954	
6	SUPPORT, CPL.		2M19654
7	SPRING	D14967	
8-10	MODIFICATION KIT CORNER PIECE, <a href="#">see picture 72</a>		2M20722
11	HANDLE, CPL., <a href="#">see picture 11</a>		2M20391
12	PROTECTING CAP (FOR AIR)	V04614	
13A	HANDWHEEL, WHITE		2M21202
13B	HANDWHEEL, ORANGE		2M19647
14	SUPPORT, LEFT		2M20692
15	SCREW DWN 562-4X20-C-ST/A2P		1338382
16	WASHER B 4,3 DIN 9021-A4		1331914
17	SUPPORT, RIGHT		2M20660
18	INCLINATION, RIGHT, <a href="#">see picture 30</a>		
19	SCREW M6x30 DIN 912-A2	1263366	
20	SPLIT WASHER B6 DIN127-STAINL.		1333941
21	B6,4 DIN 9021-A4	1292358	
22	INCLINATION, LEFT, <a href="#">see picture 30</a>		
23-24	LYING SURFACE		2M20888
23	REST SURFACE	2M20889	
24	LOCK WASHER	G12470	
25	SCREW DIN 6912-M8X20-A2	1342061	
26-27	LYING SURFACE		2M20888
26	REST SURFACE	2M20889	
27	DIVIDING WALL	2M20890	
28	MATTRESS		2M20907

Item No.	Description	Part No.	Order-Code
29	INTERMEDIATE ELEMENT		2M19537
30	TUB		2M19334
31	GASKET		2M19595
32	HOUSING NC,SC,IC, CPL., <a href="#">see picture 12 - 15</a>		
34	BED EXTENSION		2M21225
35	BORDERING,CPL., <a href="#">see picture 48</a>		
36	CABLE HARNESS		8200759
37	LABEL	2M20617	
38	ELECTRIC MAINS NC/SC VARIABLE	2M21755	
39	CABLE CLIP (2,4X92)		8712007
40	BLUE NOSE CAT without illustration		2M21420
42	MATTRESS-CLOTH / BT8000 without illustration		2M21272
43	FERRITE TUBE DI=11 LV00862	1825712	
44	FERRITE TUBE DI=8 LV00861	1825704	
45	SCREW DWN 562-4X20-C-ST/A2P	1338382	
47	CARRIAGE W. LIFTING COLUMN 120V, <a href="#">see picture 4</a>		
48	FUSE LINK DIN 41662 T6,3A		1815172
49	HEATING UNIT (120V), <a href="#">see picture 25</a>	2M21184	
50	<a href="#">See picture 72</a> and TSB # 4		

CARRIAGE WITH LIFTING COLUMN

Picture 4






Item No.	Description	Part No.	Order-Code
1-21	CARRIAGE W.LIFTING COLUMN 120V	2M21863	
1-5	HEIGHT ADJUSTMENT SC/IC 100/127V 50/60HZ		2M20244
1	ALTIT. ADJ. SC/IC 100/127V 50/60HZ		2M20296
2	PEDAL, CPL., <a href="#">see picture 5</a>		2M20290
3	SREW M5X12 DIN 912-A4	1325760	
4a	FIXING GUIDE ROLL 125 MM		2M20794
5a	GUIDE ROLL 125 MM		2M20792
4b	FIXING GUIDE ROLL 100 MM		2M20260
5b	GUIDE ROLL 100 MM		2M20250
6-21	CARRIAGE, CPL.		2M20240
6	SCREW M10 DIN 934-M A2	1335103	
7	CH. HEAD SCREW M10X40DIN912-8.8	1308564	
8	DISK DWN681-10-ST/GALV.2PIECES	1340921	
9	SCREW M6X40 DIN 912-8.8/155	1295594	
10	SCREW M6x30 DIN 912-A2	1263366	
11	WASHER B 6,4 DIN 125-A4		1330705
12	SPLIT WASHER B6 DIN127-STAINL.		1333941
13	CH. HEAD SCREW M4X8 DIN912		1263056
14	SPLIT WASHER B4 DIN 127		1331345
15	FLACHSTECKER DIN46244-A6,3 X FLAT PLUG DIN 46244-A6,3		1814079
16	CAP	2M18020	
17	BASE PLATE	2M20248	
18A	CAP, WHITE		2M21223
18B	CAP, ORANGE		2M20348
19	TUBE, BACK		2M20232
20	TUBE, FRONT		2M20234
21	SHEET METAL, CPL.	2M20584	


# **PEDAL FOR HEIGHT ADJUSTMENT**

**Picture 5**



Item No.	Description	Part No.	Order-Code
1-7	PEDAL, CPL.		2M20290
1	PEDAL		2M20218
2	PLATE, CPL.		2M20370
3	RUBBER BUFFER	M20685	
4	WASHER		G80073
5	SELF-LOCKING NUT M6 DIN 980-A2		1333062
6	B6,4 DIN 9021-A4	1292358	
7	CYL. HEAD SCREW M 6X50 DIN 912	1309129	
8	M5X12 DIN 912-A4, without illustration	1325760	


# **PLATE FOR HEIGHT ADJUSTMENT**

**Picture 6**


Item No.	Description	Part No.	Order-Code
1-7	PLATE, CPL.		2M20370
1	PLATE	2M20293	
2	MICROSWITCH		1831666
3	CABLE HARNESS	2M20213	
4	SHEET METAL	2M20365	
5	FLAT PLUG DIN 46244-A6,3		1814079
6	LOCK WASHER		1293680
7	SCREW AM 4X10 DIN84-A4/051		1329049

## FAN WHEEL


Picture 11



Item No.	Description	Part No.	Order-Code
1-3	FAN WHEEL, CPL.		2M20205
1	CONE		2M20542
2	FAN WHEEL		2M20203
3	MAGNET	2M20204	
4	TOROIDAL SEALING RING		M19241
5	SET SCREW AM 3X4 DIN 914-A2	1328255	
6	O-RING SEAL	R22364	
7	CONE (ORDER POS. 1)	2M20133	


## HOUSING NC/SC/IC

Picture 12




## HOUSING NC/SC/IC

Picture 13




HOUSING NC/SC/IC

Picture 14



HOUSING NC/SC/IC


Picture 15





Item No.	Description	Part No.	Order-Code
1-38	HOUSING NC, SC, IC, CPL.	2M21130	
	HOUSING, CPL.	2M20255	
1	SCREW DWN 562-3X12-F-A2-H	1338986	
2	SCREW DWN 562-4X20-C-ST/A2P		1338382
2a	CHEESE HEAD SCREW AM 4X16	1329855	
3	SCREW DWN 562-4X10-F-A2-H	1341227	
4	SCREW BZ 2,9X19 DIN 7971-A2		1337610
5	DIN 7982-ST2,9X9,5-F-A2-H	1341499	
6	SCREW M12x40 DIN 912-8.8/155	1264087	
7	FAN PROTECTIVE GRATING	1821997	
8	FAN SUPPORT	2M19915	
9	CONNECTOR		2M20332
10	SHEET METAL	2M20682	
11	ELBOW	2M19588	
12	LEVER	2M19613	
13	SCREW	2M19465	
14	FILTER SUPPORT	2M20289	
15	LOCK		2M20257
16	LABEL		2M20629
17	HOUSING		2M20299
18	GASKET		2M20199
19	LABEL AS OF S/N ARKA-001	2M21955	
19	LABEL UP TO S/N ARJN-00XX	2M20519	
20	HINGE		2M19730
21	LABEL	2M21768	
22	CROSS SUPPORT	8404154	
23	CABLE CLIP (2,4X92)		8712007
24	VALVE SEAT		RM08303
25	CONTROL VALVE PLATE		RM08304
26	TOROIDAL SEALING RING		D11137


Item No.	Description	Part No.	Order-Code
27	NOZZLE	R20365	
28	FILTER		8402926
29	HINGE BRACKET, CPL., <a href="#">SEE PICTURE 16</a>		
30	TAPPET, CPL., <a href="#">SEE PICTURE 17</a>		
31	SHEET METAL, CPL.	2M20284	
32	CABLE, CPL.	2M19463	
33	FAN		2M20021
34	EARTH CONDUCTOR TERMINAL		8300390
35	FERRITE TUBE DI=8 LV00861	1825704	
36	FERRITE TUBE DI=11 LV00862	1825712	
37	CAP	2M20356	
38	DIN 6799-15-X7CRNIAL 177	1341197	

**HINGE BRACKET****Picture 16**


Item No.	Description	Part No.	Order-Code
1-3	HINGE BRACKET, CPL.	2M19544	
1	HINGE BRACKET		2M19542
2	SCREW		2M19529
3	LOCK WASHER 7 DIN 6799		1309773
4	SHEET METAL	2M21392	

## TAPPET, not for Scand. Hood

Picture 17




Item No.	Description	Part No.	Order-Code
1-7	TAPPET,CPL.		2M19543
1	TAPPET	2M19532	
2	CLAMPING SLEEVE 2X8 DIN 1481	1269879	
3	LOCK WASHER 12 DIN 6799		1291467
4	THRUST WASHER 4 DIN 6799		1328999
5	GUIDE		2M19531
6	SPRING		T50137
7	TAPPET	2M19549	

**CATCH FOR SCAND. HOOD****Picture 18**

Item No.	Description	Part No.	Order-Code
1-5	SPRING-LOADED CATCH FOR SCAND. HOOD		2M20852
1	HOSE 5X2SI 60SHA NF M25779		1203606
2	BOLT		2M20853
3	LOCK WASHER 12 DIN 6799		1291467
4	GUIDE		2M19531
5	THRUST WASHER 4 DIN 6799		1328999

## SWIVEL CUPBOARD


Picture 19



Item No.	Description	Part No.	Order-Code
1-6	SWIVEL CUPBOARD, CPL.		2M20638
1	SWIVELLING DRAWER		2M20627
2	INSERT		2M20642
3	SWIVELLING PLATE		2M20628
4	TUBE		2M20639
5	CATCH		2M20632
6	SPRING		2M20645

## Mount for cupboard


Picture 20



Item No.	Description	Part No.	Order-Code
1-6	MOUNTG. DEVICES F. CUPB. INCUBATOR		2M20656
1	NUT		1267582
2	DIAPHRAGM	D15636	
3	SCREW DIN 912-M8X390+-2 8.8/155		1342096
4	DISK SPRING A25 DIN 2093		1270338
5	B8,4 DIN 9021-A4		1335367
6	SOCKET	2M20655	

## ELECTRIC MAINS NC/SC, VARIABLE HEIGHT / FIXED HEIGHT


Picture 23





## ELECTRIC MAINS NC/SC, VARIABLE HEIGHT

Picture 24




Item No.	Description	Part No.	Order-Code
1	M4 DIN 546-MS/081	1317873	
2	SPLIT WASHER B4 DIN 127		1331345
3	WASHER A 4,3 DIN 125-A2		1327542
4	SCREW DWN 562-4X20-C-ST/A2P		1338382
5	SCREW AM3X10 DIN 85-A2		1333984
6	HEXAGON NUT M3 DIN934-M A4		1334913
7	DIN 7985-M4X6-A4-H	1338862	
8	ELECTRIC MAINS 3POL X	1840339	
9	COVER	2M21724	
10	CAP without illustration	2M20356	
11	SHEET METAL	2M20697	
12	CABLE CLIP (2,4X92)		8712007
13	CROSS SUPPORT	8404154	
14	ANGLE, CPL.	2M21723	
15	CABLE	2M21764	
16	CABLE	2M21766	
17a	CABLE, AIR TEMPERATURE		8290435
17b	CABLE, SKIN TEMPERATUR (ADDITION TO 8290435)		8290445
18	SUPPORT WITH CAP LV01414	1832662	
19	CONNECTOR 6,3 LV02225	1840320	
20	CABLE	2M21739	
22	CABLE		8400808
23	SET OF CABLE GATE		2M17104
24	DISC 13 DIN 125	1330837	
25	SHEET METAL II	2M20226	
26	LABLE	8712046	
27	SCREW B 2,2X4,5 DIN 7971		1330217
27	FLAT HEAD SCREW AM 2X4 DIN84	1287664	
28	SHEET METAL	2M20224	
29	CABLE	2M20374	

Item No.	Description	Part No.	Order-Code
30	CABLE	2M20375	
31	SHR.D.PL.TUB.6,35-3,2-8404085X	1202944	
32	CONNECTOR	2M15847	
33	HEXAGON NUT M4 DIN 934-A4/051		1328956
34	CABLE	2M20385	
35	FLAT PLUG DIN 46244-A6,3		1814079
36	CAP	M27256	
37	SCREW DWN 561-4X20-C-A2-H		1338072
38	SCREW AM3X10 DIN 85-A2		1333984
39	SPLIT WASHER B3 DIN127-X12CRNI		1309749
40	ADH. TAPE 12X1,6 PVC-FOAM SW X		1202421
41	LABEL	2M20631	
43	MAINS CABLE, UL		2M20453
44	MAINS CABLE BUSHING		2M18642
45	SICHERUNGSEINS.DIN41662 T6,3A FUSE LINK DIN 41662 T6,3A		1815172

## AGGREGATE

Picture 25




Item No.	Description	Part No.	Order-Code
1-49	AGGREGATE, CPL. with humidifier, heating cartridges, fuses, fan wheel		2M20615
1-49	REP.EXCH.AGGREG.HEAT. (2M20615)	2M20588	
1	GUIDE	2M20346	
2	CABLE HARNESS	2M20357	
3	CROSS SUPPORT	8404154	
4	M4 DIN 546-MS/081	1317873	
5	DISSIPATOR	2M19615	
6	SCREW AM4X6 DIN 84-A4		1330829
7	SPLIT WASHER B4 DIN 127		1331345
8	CABLE CLIP (2,4X92)		8712007
9	BOLT	2M20347	
10	WASHER B 4,3 DIN 9021-A4		1331914
11	HEXAGON NUT M4 DIN 934-A4/051		1328956
12	TEMPERATURE FUSE 315 DEGREE		1823809
13	SILICONE HOSE 6X2,5 NF M29907		1197851
14	SCREW AM 4X10 DIN84-A4/051	1329049	
15	DISC B 5,3 DIN 9021-A4		1329472
16	SCREW AM 4X12 DIN 84-A2		1287710
17	SCREW NUT	2M20297	
18	PCB SHEET METAL, CPL.	2M20325	
19	CABLE HARNESS	2M20369	
20	GASKET		2M19637
21	SHEET METAL, CPL.	2M20322	
22	HEXAGON NUT M6 DIN 934-M A4	1267515	
23	A6 DIN 137-STAINLESS STEEL		1332015
24	MOTOR,CPL., <a href="#">SEE PICTURE 26</a>		
25	SPLIT WASHER B3 DIN127-X12CRNI		1309749
26	HEXAGON NUT M3 DIN934-M A4		1334913
27	WASHER		1330233

Item No.	Description	Part No.	Order-Code
28	SCREW DIN 84-AM3X20-4.8/155	1303325	
29	PACKING RING		2M19638
30	PLATE	2M20023	
31	COIL		2M18839
32	PROTECTING CAP	D13337	
33	HUMIDIFIER, <a href="#">SEE PICTURE 28</a>		
34	ADHESIVE LABEL, without illustration		8303296
37	FAN WHEEL, CPL., <a href="#">SEE PICTURE 11</a>		
38	HEATING CARTRIDGE (120/127V) FOR HUMIDIFIER, <a href="#">SEE PICTURE 28</a>		
39	HEATING CARTRIDGE 120/127 V		2M20331
40 - 42	Fuse T 6,3 A DIN 41662		1815172
40 - 42	Fuse T 1A DIN41662		1815083
40 - 42	Fuse T 80 mA DIN 41662		1814974
40 - 42	Fuse T 200 mA DIN 41662		18 15 016
40 - 42	Fuse T 500 mA DIN 41662		18 15 059
40 - 42	Fuse T 2,5 A DIN 41662		18 15 121
40 - 42	Fuse T 1,6 A DIN 41662		18 15 105
40 - 42	Fuse T 160 mA DIN 41662		18 15 008
46	CONNECTION, without illustration		2M20464
50	CABLE HARNESS (PLUG-IN UNIT - AGGREGATE) without illustration		8200759

## PCB SHEET METAL (PART OF AGGREGATE)

Picture 26




Item No.	Description	Part No.	Order-Code
1-21	PCB SHEET METAL, CPL.	2M20325	
1	RELAY SEMICONDUCTOR		1831399
2	SCREW AM 4X10 DIN84-A4/051	1329049	
3	SPLIT WASHER B4 DIN 127		1331345
4	WASHER A 4,3 DIN 125-A2		1327542
5	EARTH CONDUCTOR TERMINAL		8300390
6	SCREW M3X6 DIN 84-A4		1330810
7	SPLIT WASHER B3 DIN127-X12CRNI		1309749
8	SHEET METAL, CPL.	2M20355	
9	TRANSFORMER		8290426
10	SCREW DIN 84-AM4X50-A4		1287753
11	CABLE CLIP (2,4X92)		8712007
12	HEADER PROFILE 1-2CH15956 SW		1191640
13	CABLE CLAMP	2M20380	
14	SCREW		1329987
15	WASHER		1330233
16	PCB SHEET METAL, CPL.	2M20324	
17	PCB AGGREGATE		8290511
17a	RELAY 12V/8A 1XUM, (K1/K5 ON AGGREGATE PCB)		1830732
17b	RELAY 24V/6A/250V, (K2/K5 ON AGGREGATE PCB)		1821474
17c	RELAY, (K4 ON AGGREGATE PCB)		8301198
18	CABLE BUSHING	CH17949	
19	PCB HL-RELAY		8306151
20	B5 DIN 137-STAINLESS STEEL	1335316	
21	B4 DIN 137-STAINLESS STEEL		1335308
22	DISTANCE BOLT *M3*, without illustration	8304504	



MOTOR

Picture 27




Item No.	Description	Part No.	Order-Code
1-22	MOTOR,CPL. (ORDER IN ADDITION CONE 2M20542 FOR MOTOR REPLACEMENT)		2M20350
1	LIP SEAL	2M20582	
2	V-SHAFT SEAL	2M11477	
3	CONNECTION	2M19689	
4	SHEET METAL	2M19676	
5	MOTOR . (ORDER IN ADDITION CONE 2M20542 FOR MOTOR REPLACEMENT)		2M20379
6	SHEET METAL	2M19677	
7	MOTOR SUPPORT	2M20616	
8	GASKET		2M19591
9	RUBBER SUSPENSION		2M18147
10	NUT	2M19675	
11	SPRING		E20215
12	CABLE CLIP (2,4X92)		8712007
13	NUT BM3 DIN 439-A4/051		1334743
14	COUNTER SUNK AM4X16DIN963-A2	1329332	
15	SCREW DWN 561-4X20-C-A2-H		1338072
16	SET SCREW M3X8 DIN 551-A2	1324926	
17	SCREW M3X25 DIN 933-A2	1341715	
18	HEXAGON NUT M3 DIN934-M A4		1334913
19	SHR. D. PL. TUB. 6,35-3,2-8404085X	1202944	
20	HOSE 6732566, 3X0,75	1201778	
21	CABLE CLAMP	2M20380	
22	CABLE	2M20298	
23	BOTTLE OF SPECIAL OIL, without illustration		2M07839
24	SCREW AM6X10 DIN 963-A2	1334646	
25	BOLT	2M19678	
26	BOLT	2M19679	
27	MOTOR SUPPORT, for exchange please order item no. 1-22		

Item No.	Description	Part No.	Order-Code
28	CROSS SUPPORT	8404154	
29	SCREW AM4X6 DIN 84-A4		1330829
30	WASHER A 4,3 DIN 125-A2		1327542
31	SPLIT WASHER B4 DIN 127		1331345

## HUMIDIFIER


Picture 28



Item No.	Description	Part No.	Order-Code
1-14, 16,18	HUMIDIFIER, CPL., 120/127V		2M20619
1	CONNECTION		2M20292
2	O-RING SEAL		2M08777
3	XFL 19POL STGE CRIMP LV01587	1831631	
4	HOSE CLIP	T12881	
5	SCREW AM 4X12 DIN 84-A2		1287710
6	SPLIT WASHER B4 DIN 127		1331345
7	WASHER A 4,3 DIN 125-A2		1327542
8	TAPPET	2M20359	
9	CABLE	2M20383	
10	COOKER		2M20291
11	TEMPERATURE LIMITER		2M20382
12	TEMPERATURE SENSOR, CPL.		2M20381
13	CHEESE HEAD SCREW AM 4X16	1329855	
14	SCREW M4X4 DIN 916-A4	1306626	
16	HEATING CARTRIDGE (120/127V)		2M20328
18	SCREW AM4X6 DIN 84-A4		1330829

## HANDLE


Picture 29



Item No.	Description	Part No.	Order-Code
1-5	HANDLE, CPL.		2M20391
1	HANDLE	2M20698	
2	END PIECE		2M20696
3	TUBE		2M13388
4	SCREW M 8X55 DIN 912-A2		1329693
5	CAP SCREW M 8X75 DIN 912	1326473	
3,5-9	HANDLE,CPL.		2M20391
6	HANDLE	2M20392	
7	B8,4 DIN 9021-A4		1335367
8	SCREW M8X60 DIN 912-A2	1341286	
9	CLOSING PLUG		G60036


## INCLINATION, RIGHT, LEFT

Picture 30




Item No.	Description	Part No.	Order-Code
1-19	INCLINATION, RIGHT		2M19712
1-19	INCLINATION, LEFT		2M19650
1-10	COUPLING	2M19713	
1	LOCKING PISTON	2M19643	
2	BEARING	2M19649	
3	4X22 DIN 1481-X12CRNI177	1329758	
4	CHAIN WHEEL	2M19661	
5	CLAMPING SLEEVE 3X12 DIN 1481	1269968	
6	SPRING		2M19648
7	AXLE	2M19644	
8	WASHER		CH09520
9	SEALING RING 22X1 DIN 472	1297082	
10	PIN	2M19651	
11	BASE PLATE	2M19646	
12	CHAIN		2M19660
13	CHAIN WHEEL DISK	2M19662	
14	SCREW	2M19652	
15	SCREW	2M19653	
16	PLATE	2M19708	
17	607 DIN 625	1330128	
18	RECEIVER, CPL.		2M19657
19	SCREW AM 2X8 DIN 84-A4	1327984	



**MODIFICATION SET H2O / INC 8000 SC TO IC****Picture 31**

## MODIFICATION SET H2O / INC 8000 SC TO IC


Picture 32



Item No.	Description	Part No.	Order-Code
1-27	MODI. SET H2O / INC. 8000 SC TO IC		2M20820
1	SCREW DWN 561-4X20-C-A2-H, without illustration		1338072
2	BUFFER		6800478
3	BUTTON		2M19693
4	INFUSION BOTTLE DIN 58363 (NOT FOR US VERSION)	1340697	
5	SCREW M4X12 DIN 427-A2		1332007
6-21	FLAP 8000 IC, CPL.	2M20455	
6	RETAINING NUT M5 DIN 980-A2	1335618	
7	HEXAGON NUT M4 DIN 985	1298720	
8	DISC B 5,3 DIN 9021-A4		1329472
9	WASHER B 4,3 DIN 9021-A4		1331914
10	14,5/20 DIN 12871 SILICONE, without illustration	1291920	
11	BOLT	2M19533	
12	GASKET (USA/CAN)		2M20635
13	KNURLED NUT		2M06004
14	BUTTON		2M12044
15	TRACTION ROPE		2M12045
16	SPRING-LOADED CATCH		2M19271
17	SUPPORT	2M19806	
19	WATER TANK		2M19561
20	CONNECTION		2M20464
21	FLAP	2M20463	
22-27	SHEET METAL		2M20461
22	CAP		2M20459
23	FLAT PLUG DIN 46244-A6,3		1814079
24	SPLIT WASHER B4 DIN 127		1331345
25	HEXAGON NUT M4 DIN 934-A4/051		1328956
26	ADH. TAPE 12X1,6 PVC-FOAM SW X		1202421
27	SHEET METAL		2M20549

WATER CONTAINER SET


Picture 33



Item No.	Description	Part No.	Order-Code
1-11	WATER CONTAINER SET		2M21747
1-6	WATER CONTAINER, CPL.	2M20272	
1,7-11	SHEET METAL, CPL.	2M21642	
1	SCREW DWN 561-4X20-C-A2-H		1338072
2-5	WATER CONTAINER,SET OF SPARES		2M20644
2	CONNECTION		2M20670
3	WATER TANK		2M20210
4	COVER		2M20211
5	LABEL	2M20562	
6	CLAMP		2M20547
7	SPLIT WASHER B4 DIN 127		1331345
8	HEXAGON NUT M4 DIN 934-A4/051		1328956
9	FLAT PLUG DIN 46244-A6,3		1814079
10	SHEET METAL	2M20224	
11	ADH. TAPE 12X1,6 PVC-FOAM SW		1202421

STANDARD HOOD W/O DOUBLEWALLS (OLD VERSION)


Picture 34



Item No.	Description	Part No.	Order-Code
1-18	HOOD INCUBATOR 8000 SC, CPL., <a href="#">see picture 42</a>	2M22055	
1A	HOOD 8000 SC, <a href="#">see picture 42</a>	2M22038	
1B	SET COMBINATION STANDARD HOOD IC, SC, NC		2M22153
2	FLAP, CPL., <a href="#">see picture 35</a>		
3	CATCH, CPL.		2M20358
4	SCREW DWN 562-3X16-F-A2-H	1342630	
5	SET OF SWIVELLING WINDOWS, <a href="#">see picture 36</a>		
6	SWIVELLING WINDOW, LEFT, <a href="#">see picture 36</a>		
7	SWIVELLING WINDOW, RIGHT, <a href="#">see picture 36</a>		
8	CABLE CLIP		2M20376
9	SCREW DWN 561-3,5X16-C-A2	1339311	
10	CAP NUT		2M19484
11	SCREW DWN 562-3X12-F-A2-H	1338986	
12	HOSE THROUGH-BEARING		2M19511
13	STIFFENING		2M20214
14	ANGLE, <a href="#">see picture 72</a>		
15	SCREW DWN 562-5X12-F-A2-H	1342649	
16	STIFFENING	2M20568	
17	SCREW DWN 562-3X10-F-A2-H	1342622	

FLAP FOR STANDARD HOOD

Picture 35






Item No.	Description	Part No.	Order-Code
1-13	FLAP, CPL.		2M19525
1	SCREW DWN 562-3X12-F-A2-H	1338986	
2-5	CATCH, CPL.		2M19899
2	CATCH		2M19482
3	HOUSING		2M19483
4	SPRING		8400496
5	PLATE		2M19481
6	CAP NUT		2M19484
	PLATE FOR CATCH 2M19899		2M19472
7	CATCH		2M19516
8	HINGE		2M19729
9	INTERIOR WALL		2M19512
10	FLAP		2M20733
11	SCREW DWN 562-4X14-F-A2-H		1342568
12	PLATE	2M20731	
13	HINGE		2M20730

## SET OF SWIVELLING WINDOWS


Picture 36



Item No.	Description	Part No.	Order-Code
3,5-18	SET OF SWIVELLING WINDOWS, COMPLETE, LEFT AND RIGHT		2M19550
1-3	FLAP, LEFT (WITH HOSE BUSHING)		2M20552
2-4	FLAP, RIGHT (WITH HOSE BUSHING)		2M20545
1	FLAP, LEFT	2M20433	
2	TUBE BUSHING		2M20434
3	SCREW M6X6 DIN551-PA 6.6 NF		1342266
4	FLAP, RIGHT	2M20546	
5	CATCH		2M19519
6	SCREW BZ 2,9X19 DIN 7971-A2		1337610
7	SCREW BZ 2,9X16 DIN7971-A2		1336630
8	BASE PLATE		2M19521
9	SUPPORT, TOP LEFT		2M19473
10	FLAP, LEFT		2M19477
11	FLAP, RIGHT If item 10 or 11 are ordered singly, item 3, screw 1342266, must also be ordered		2M19478
12	SUPPORT, TOP RIGHT		2M19475
13	BUFFER		2M19524
14	SUPPORT, BOTTOM RIGHT		2M19476
15	SUPPORT, BOTTOM LEFT		2M19474
16	GASKET, LEFT		2M19469
17	GASKET, RIGHT		2M19470
18	SPRING		2M20537
19	SCREW	1341642	
20	GASKET (SWIVELLING WINDOW, LEFT)	2M20953	

STANDARD HOOD W/O DOUBLEWALLS (OLD VERSION)


Picture 37



Item No.	Description	Part No.	Order-Code
1-16	HOOD INC. 8000 SC, CPL., without double walls, <a href="#">see picture 42</a>	2M20220	
1A	PLEXIGLASS HOOD INC. 8000 SC		2M20196
1B	Combination Standard Hood IC, SC, NC		2M22153
2	FLAP, CPL., <a href="#">see picture 35</a>		2M19525
3	SET OF SWIVELLING WINDOWS, <a href="#">see picture 36</a>		2M19550
4	CATCH, RIGHT		2M19506
5	SCHRAUBE BZ2,9X16 DIN7971-A2		1336630
6	RASTE, LINKS		2M19505
7	SCREW BZ 2,9X13 DIN 7971-A2	1334948	
8	KABELOESE		2M20376
9	CAP NUT		2M19484
10	SCREW DWN 561-3,5X16-C-A2	1339311	
11	SCREW DWN 562-3X12-F-A2-H	1338986	
12	STIFFENING		2M20214
13	HOSE THROUGH-BEARING		2M19511
14	ANGLE, <a href="#">see picture 72</a>		
15	SCREW DWN 562-5X12-F-A2-H	1342649	
17	ELBOW, <a href="#">see picture 72</a>		

SCAND. HOOD W/O DOUBLEWALLS (OLD VERSION)


Picture 38



Item No.	Description	Part No.	Order-Code
1-21	HOOD INCU.8000 SC SCAND. HOOD WITHOUT DOUBLE WALLS, <a href="#">see picture 43</a>		2M20830
1	HOOD SC, SCAND., <a href="#">see picture 43</a>	2M20831	
2	FLAP, BIG (CPL.), <a href="#">see picture 39</a>		
3	SWIVELLING WINDOW, RIGHT, <a href="#">see picture 36</a>		
4	CATCH, CPL.		2M20358
5	SCREW DWN 562-3X16-F-A2-H		1342630
6	SWIVELLING WINDOW, LEFT, <a href="#">see picture 36</a>		
7	SCREW DWN 562-3X12-F-A2-H	1338986	
8	FLAP, SMALL (CPL.), <a href="#">see picture 40</a>		
9	SET OF SWIVELLING WINDOWS, <a href="#">see picture 36</a>		
10	CABLE CLIP		2M20376
11	CAP NUT		2M19484
12	SCREW DWN 561-3,5X16-C-A2	1339311	
13	HOSE THROUGH-BEARING		2M19511
14	STIFFENING		2M20214
15	ANGLE, <a href="#">see picture 72</a>		
16	SCREW DWN 562-5X12-F-A2-H	1342649	
17	STIFFENING	2M20568	
18	SCREW DWN 562-3X10-F-A2-H	1342622	
19	SPRING-LOADED CATCH, without illustration	2M20852	
20	SET SCREW M 5X16 DIN 427	1324004	


FLAP BIG, FOR SCAND. HOOD

Picture 39






Item No.	Description	Part No.	Order-Code
1-14	FLAP, BIG (CPL.)		2M20607
1	BUFFER		2M19524
2	FLAP, BIG		2M20604
3	CATCH		2M19516
4	SCREW DWN 562-3X12-F-A2-H	1338986	
5	CAP NUT		2M19484
6-9	CATCH,CPL.		2M19899
6	PLATE		2M19481
7	SPRING		8400496
8	HOUSING		2M19483
9	CATCH		2M19482
10	HINGE		2M19729
11	SET OF SWIVELLING WINDOWS, <a href="#">SEE PICTURE 36</a>		
12	INTERIOR WALL, BIG	2M20605	
13	HINGE	2M20606	
14	SCREW DWN 562-3X16-F-A2-H	1342630	

**FLAP SMALL, FOR SCAND. HOOD****Picture 40**

**FLAP SMALL, FOR SCAND. HOOD**


**Picture 41**



Item No.	Description	Part No.	Order-Code
1-9	FLAP, SMALL (CPL.)		2M20612
7,10-15	FLAP, SMALL (CPL.)		2M20612
1	SCREW DWN 562-4X14-F-A2-H	1342568	
2	FLAP, SMALL		2M21471
3	CATCH		2M19516
4	SCREW DWN 562-3X16-F-A2-H	1342630	
5	INTERIOR WALL, SMALL	2M21472	
6	PLATE	2M20731	
7	HINGE		2M19729
8	SCREW DWN 562-3X12-F-A2-H	1338986	
9	HINGE	2M21457	
10	FLAP, SMALL	2M20608	
11	BOLT	2M20613	
12	BZ 2,9X13 DIN 7971-A2		1334948
13	INTERIOR WALL, SMALL	2M20609	
14	O-RING SEAL		R50543
15	AXLE		2M19848

STANDARD HOOD DOUBLEWALLS (NEW VERSION)


Picture 42



Item No.	Description	Part No.	Order-Code
1-24,27	HOOD INCUBATOR 8000 SC DOUBLEWALL, COMPLETE		2M21097
1	SCREW DWN 562-3X10-F-A2-H	1342622	
2	SCREW DWN 562-3X16-F-A2-H	1342630	
3	SCREW DWN 562-3X12-F-A2-H	1338986	
4	SCREW DWN 561-3,5X16-C-A2	1339311	
5	SCREW DWN 562-5X12-F-A2-H	1342649	
6	STIFFENING	2M20568	
7	BOLT FOR REMOVEABLE INNER WALLS		2M21448
8	RECEIVER FOR REMOVEABLE INNER WALLS		2M21449
9	CABLE CLIP		2M20376
10	HOSE THROUGH-BEARING		2M19511
11	COMBINATION PLEXIGALS HOOD IC/SC/NC		2M22153
12	ANGLE, <a href="#">see picture 72</a>		
13	STIFFENING		2M20214
14	CAP NUT		2M19484
15	CATCH		2M19724
16	INTERIOR WALL		2M19512
17	LATERAL WALL		2M20508
19	CATCH, CPL.		2M20358
20	FLAP, CPL., <a href="#">see picture 35</a>		
21	SET OF SWIVELLING WINDOWS, <a href="#">see picture 36</a>		
22	SWIVELLING WINDOW, LEFT, <a href="#">see picture 36</a>		
23	SWIVELLING WINDOW, RIGHT, <a href="#">see picture 36</a>		
25	CATCH, LEFT		2M19505
26	CATCH, RIGHT		2M19506
27	LABEL, without illustration	2M21596	

## SCAND. HOOD DOUBLEWALLS (NEW VERSION)

Picture 43




Item No.	Description	Part No.	Order-Code
1-29	HOOD INCU.8000 SC WITH DOUBLE WALLS, COMPLETE		2M21620
1	SCREW DWN 562-5X12-F-A2-H	1342649	
2	SET SCREW M 5X16 DIN 427	1324004	
3	SCREW DWN 562-3X16-F-A2-H	1342630	
4	SCREW DWN 562-3X12-F-A2-H	1338986	
5	SCREW DWN 561-3,5X16-C-A2	1339311	
6	SCREW DWN 562-3X10-F-A2-H	1342622	
8	STIFFENING		2M20568
9	HOOD SC SCAND.W.DOUBLE WALL		2M21619
10	CAP NUT		2M19484
11	INTERIOR WALL		2M19512
12	CABLE CLIP		2M20376
13	HOSE THROUGH-BEARING		2M19511
14	STIFFENING		2M20214
15	ANGLE, <a href="#">see picture 72</a>		
16	SIDE WALL		2M21727
17	CATCH		2M19724
18	RECEIVER FOR REMOVEABLE INNER WALLS		2M21449
19	BOLT FOR REMOVEABLE INNER WALLS		2M21448
20	FLAP, BIG (CPL.), <a href="#">see picture 39</a>		
21	SET OF SWIVELLING WINDOWS, <a href="#">see picture 36</a>		
22	FLAP, SMALL (CPL.), <a href="#">see picture 40</a>		
23	SWIVELLING WINDOW, LEFT, <a href="#">see picture 36</a>		
24	SWIVELLING WINDOW, RIGHT, <a href="#">see picture 36</a>		
25	SPRING-LOADED CATCH, without illustration		2M20852
26	CATCH, CPL.		2M20358



SCAND. HOOD W/O DOUBLEWALLS (OLD VERSION)


Picture 46



Item No.	Description	Part No.	Order-Code
1-21	HOOD INCU.8000 SC SCAN. HOOD, <a href="#">see picture 43</a>	2M22110	
1	HOOD SC, SCAN, <a href="#">see picture 43</a>	2M22109	
2	FLAP, BIG (CPL.), <a href="#">see picture 39</a>		
3	SWIVELLING WINDOW, RIGHT, <a href="#">see picture 36</a>		
4	CATCH, CPL.		2M20358
5	SCREW DWN 562-3X16-F-A2-H	1342630	
6	SWIVELLING WINDOW, LEFT, <a href="#">see picture 36</a>		
7	SCREW DWN 562-3X12-F-A2-H		1338986
8	FLAP, SMALL (CPL.), <a href="#">see picture 39</a>		
9	SET OF SWIVELLING WINDOWS, <a href="#">see picture 36</a>		
10	CABLE CLIP		2M20376
11	CAP NUT		2M19484
12	SCREW DWN 561-3,5X16-C-A2	1339311	
13	HOSE THROUGH-BEARING		2M19511
14	STIFFENING		2M20214
15	ANGLE, <a href="#">see picture 72</a>		
16	SCREW DWN 562-5X12-F-A2-H	1342649	
17	STIFFENING		2M20568
18	SCREW DWN 562-3X10-F-A2-H	1342622	
19	SPRING-LOADED CATCH, without illustration		2M20852
20	SET SCREW M 5X16 DIN 427	1324004	

## TEMPERATURE SENSOR


Picture 47



Item No.	Description	Part No.	Order-Code
1-6	TEMPERATURE-SENSOR INC.8000 SC		2M20277
1	SCREW BZ 2,9X16 DIN7971-A2		1336630
2	COVER	2M20334	
3	SUPPORT	2M20335	
4	TEMPERATURE SENSOR	2M20282	
5	SHRINKD.TB.1,2-0,6-NF 6803472X	1203819	
6	CABLE	2M20386	

## BORDERING


Picture 48



Item No.	Description	Part No.	Order-Code
1	PROFILE		2M20195
2	BORDERING, WHITE	2M21201	
3	WASHER A 4,3 DIN 125-A2		1327542
4	LOCKING PIN		2M20351
5	SCREW DWN 562-4X20-C-ST/A2P		1338382
6	SHEET METAL	2M20338	
7	FLAT PLUG DIN 46244-A6,3		1814079
8	SPLIT WASHER B4 DIN 127		1331345
9	HEXAGON NUT M4 DIN 934-A4/051		1328956
10	ADH. TAPE 12X1,6 PVC-FOAM SW X		1202421
11	LABEL, without illustration		2M20564
12	BUFFER		2M21887


## PLUG-IN UNIT 8000 NC/SC

Picture 49




PLUG-IN UNIT 8000 NC/SC

Picture 50




**PLUG-IN UNIT 8000 NC/SC**

**Picture 51**




30.01.1996


**PLUG-IN UNIT 8000 NC/SC****Picture 52**



PLUG-IN UNIT 8000 NC/SC


Picture 53



**PLUG-IN UNIT 8000 NC/SC****Picture 54**


PLUG-IN UNIT 8000 NC/SC

Picture 55



Item No.	Description	Part No.	Order-Code
	PLUG-IN UNIT 8000 SC with Skin Temp. Control, without EPROM, for SW 10.n only, without central alarm LED		8290441
	SERVICE EXCHANGE PLUG-IN UNIT (8290441)	8290603	
1	NUT BM3 DIN 439-A4/051		1334743
2	HEXAGON NUT M3 DIN934-M A4		1334913
3	SPLIT WASHER B3 DIN127-X12CRNI		1309749
5	ADH.TAPE 19X2-TESAMOLL 4700 WS		1209582
5	ADH.TP.19 TESABD651SW,REEL 25M	1186965	
6	BUTTON		8200619
7	SIDE PART	8290408	
8	SIDE PANEL	8200893	
9	BRACE	8200895	
10	BRACE	8200896	
11	SOCKET	8200897	
12	FIXING BAR		8200899
13	PLATE	8200916	
14	BELLOW		8200927
15	SLEEVE	8200116	
16	FIBRE RING		D04766
17	INSULATING BUSH		6804141
18	T-K-RAM, <a href="#">see picture 58</a>		
19	INSULATING WASHER	1814303	
20	FLAT PLUG DIN 46244-A6,3		1814079
21	COUNTERSUNK SCREW AM 3X4DIN963	1308459	
22	SCREW		1287419
23	COUNTERS. SCREW AM3X6DIN963-A4		1287400
24	SCREW		1329987
25	CYLIND. HEAD SCREW AM4X20 DIN84	1328980	
26	SCREW AM3X16 DIN 963-A2		1334557
27	WASHER A 4,3 DIN 125-A2		1327542


Item No.	Description	Part No.	Order-Code
28	SERR. LOCK WASHER A3,2 DIN6798	1270907	
29	CH.EAD SCR. AM3XDIN84-A4/051		1330799
30	PCB DISPLAY AIR, <a href="#">see picture 56</a>		
31	PCB MOTHERBOARD		8290471
32	D PIA 6821 1MHZ DIL 40 LV00683		1823361
33	CABLE HARNESS	8200759	
34	PCB SWITCH		8200880
35	PCB POWER PACK, <a href="#">see picture 57</a>		
36	PCB CPU INCUBATOR, <a href="#">see picture 58</a>		
37	CABLE HARNESS,VOLTAGE REGULAT.		8306259
39	PCB DISPLAY HUMIDITY, <a href="#">see picture 60</a>		
40b	PCB DISPLAY for Skin. Temp. Module with SW 10.n only, <a href="#">see picture 61</a>		
45	CAP		8200615
46	NREGF 5V TO220 LV01257		1830724
47	SHR. D. PLAST. TUB. 3,2-1,6X0,5	1200186	
48	RELAY BIPO 24/2 (part of Motherboard PCB, without illustration)		1812653
49A	KEY PAD, CPL., SC, with Skin Temp. Control, for SW 10.n only, without central alarm LED (replaced by P/N 2M 21 691)		8290415
49B	KEY PAD, CPL., SC, with Skin Temp. Control, for SW 10.n only, with central alarm LED		2M21691
51	CABLE HARNESS,VOLTAGE REGULAT.		8306259
52, 53	PCB ANALOG, <a href="#">see picture 59</a>		
54	CABLE, SKIN TEMPERATURE		8290445
55	PCB CONTROLLER without illustration, <a href="#">see picture 62</a>		

**PCB DISPLAY AIR****Picture 56**

Item No.	Description	Part No.	Order-Code
	PCB DISPLAY AIR		8290411
D1	D PROG.KEY LV00376		1820877

## PCB POWER PACK


Picture 57



Item No.	Description	Part No.	Order-Code
1	PCB POWER PACK		8290461
2	ACOUSTIC SIGNAL GENERATOR		8301708
3	CABLE CLIP (3,6X140)	8300358	
4	NC STORAGE BATTERY 9V		8301856
5	ADDITIONAL FAN (not necessary)		See TSB # 17

## PCB CPU INCUBATOR

Picture 58




Item No.	Description	Part No.	Order-Code
	PCB CPU INCUBATOR		8290571
D1	INTEGRATED CIRCUIT MC 6802	1809156	
D6	MEMORY DECODER, CPL.		8304920
D8	D-RAM 16K+BATT+CLOCK (MK48T02B)		1828142
D12	SOFTWARE 8000 SC/NC 10.05		2M22326
D13	TIMER 6840		1813900
D17	MONITORING LOGIC,CPL.	8304981	
X1	XFL SHORT CIRCUIT SOCKET 2POL.		1814508




## PCB ANALOG

Picture 59




Item No.	Description	Part No.	Order-Code
PCB Analog for Incubator 8000 SC without ThermoMonitoring and with SW 10.n:			
	PCB ANALOG with Skin Temp. Control, for all version with SW 20.04 (replaces P/N 8200920, see TSB # 20)		2M22381
	PCB ANALOG W/O THERMOMONITORING, as of SW version 10.05	2M22403	
	PCB ANALOG W/O THERMOMONITORING, for all SW versions 10.n	2M22402	
	PCB ANALOG W/O THERMOMONITORING, for all SW versions 10.n	2M22401	
	PCB ANALOG W/O THERMOMONITORING, for all SW versions 10.n	2M22400	
K1,K2	RELAY 5V/0,5A (K1/K2 of Analog PCB P/N 2M22400 and P/N 2M22401)		1813714
PCB Analog for Incubator 8000 IC <b>without</b> ThermoMonitoring and with SW 11.n:			
	PCB ANALOG, THERMOMONITORING, for all versions with SW 11.02 (replaces P/N 8290680, see TSB # 20)		2M22382

Item No.	Description	Part No.	Order-Code
	PCB ANALOG WITH THERMOMONITORING, AS OF SW 11.02	2M22407	
	PCB ANALOG WITH THERMOMONITORING, AS OF SW 11.02	2M22406	
	PCB ANALOG WITH THERMOMONITORING, FOR ALL SW VERSIONS 11.N	2M22405	
	PCB ANALOG WITH THERMOMONITORING, FOR ALL SW VERSIONS 11.N	2M22404	

**PCB DISPLAY HUMIDITY****Picture 60**


Item No.	Description	Part No.	Order-Code
1	PCB DISPLAY HUMIDITY		8290451
2	D PROG.KEY		1820877

**PCB DISPLAY****Picture 61**

Item No.	Description	Part No.	Order-Code
1	PCB DISPLAY		8290691
2	D PROG.KEY		1820877

## INSTALLATION SET RS 232


Picture 62



Item No.	Description	Part No.	Order-Code
	INSTALLATION SET IC8000 RS 232		8290607
	PCB CONTROLLER	8290581	
D17/D19	SOFTWARE 1.00 BABYLINK	8290599	
	CABLE HARNESS RS 232 INCUBATOR without illustration	8290568	
	BABYLINK CABLE BL.8000 without illustration		8306488

## CABLE HARNESS SKIN TEMPERATURE MEASUREMENT


Picture 63



Item No.	Description	Part No.	Order-Code
1-11	CABLE HARNESS		2M21738
1	SHEET METAL	2M21694	
2	AM6X30 DIN 963-A2	1334689	
3	CABLE HARNESS SENSOR CONNECTOR YELLOW.	2M21678	
4	SCREW M3X6 DIN 912-A2	1331930	
5	XFL 7POL BUGE LV00271		1817604
6	SOCKET BOARD COD	1813536	
7	LABEL		2M21759
8	SCREW DWN 562-4X10-F-A2-H	1341227	
9	CROSS SUPPORT	8404154	
10	FERRITE TUBE DI=8 LV00861	1825704	
11	CABLE CLIP (2,4X92)		8712007

# SKINTEMPERATURE SENSORS

Picture 66






Item No.	Description	Part No.	Order-Code
For Incubator 8000 SC/NC with connection for the skin temperature sensor on the left side of the blue housing:			
1	SENSOR-SKIN TEMP-YELLOW (5 PACK), NEW VERSION, FOR CENTRAL SKIN TEMPERATURE		MX11000
	HYDRO TEMP PROBE COVERS (100 pack) without illustration		DHNICU-00
For Incubator 8000 SC/NC with connection for the skin temperature sensor on the back side of housing:			
6	SKIN-TEMP. SENSOR, (NON DISPOSABLE), CAN BE USED WITHOUT ADAPTER 2M20736 NOTE: DISCONTINUING PART	8200757	
7	ADHESIVE RING, 250 PIECES, FOR P/N 8200757		6850021
9	INTERFACE ADAPTER CABLE / SENSOR MALLINC.		2M20736
2	TEMP SENSOR, SINGLE USE, BLUE, FOR CENTRAL TEMPERATURE (WITH ADHESIVE PAD), TO BE USED WITH ADAPTER CABLE P/N 2M20736		D5030170
	HYDRO TEMP PROBE (BLUE ANGLED CONNECTOR), CAN BE USED WITHOUT ADAPTER 2M20736		DHNICU-32
	HYDRO TEMP PROBE COVERS (100 pack) without illustration		DHNICU-00

## MODIFICATION SET LANGUAGE


Picture 67



Item No.	Description	Part No.	Order-Code
1-6,9-12	MODIFICATION SET INCU. SC USA	2M21948	
1-2	OPERATING INSTR.(INCUB.8000 SC	2M21906	
1	FLAP (INCUBATOR 8000 SC)	2M21968	
2	RUBBER BUFFER	2M21865	
3	LABEL, ENGL./USA(OPERAT.INST.SC	2M21878	
11	LABEL	2M20560	

## VENTILATION HOSE HOLDER


Picture 68



Item No.	Description	Part No.	Order-Code
1-9	VENTILATION HOSE HOLDER, new		8411075
1,5,10,11	VENTILATION HOSE HOLDER		8411075
1	SCREW		2M15332
2	CONE	8411018	
3	FLEXIBLE BRACKET	8410065	
4	CONNECTION	8411017	
5	SCREW B 2,9X19 DIN7973	1328654	
6	CATCH	8411077	
7	CLAMP		8411076
8	DISC		1335359
9	SCREW BZ 2,9X13 DIN 7971-A2		1334948
10,11	SET OF TUBE CLIPS 2		8406138
10	TUBE CLIP 2	8404761	
11	HOSE SUPPORT	2M19670	

## MODIFICATION SET USA


Picture 69



Item No.	Description	Part No.	Order-Code
1-10	MODIFICATION SET INCUBATOR 8000 IC USA	2M21984	
1-2	DEPOSIT TRAY, COMPL.		2M20320
1	DEPOSIT PLATE	2M20695	
2	RUBBER BUFFER	2M20704	
3	GASKET		2M20635
4	LABEL	2M20631	
5	FLAT HEAD SCREW AM 2X4 DIN84		1287664
6-10	SET OF LABELS USA		2M20580
6	LABEL	2M20557	
7	LABEL	2M20558	
8	LABEL	2M20559	
9	LABEL	2M20556	
10	LABEL	2M20561	

## MODIFICATION KIT INCU. 8000 SC USA


Picture 70



Item No.	Description	Part No.	Order-Code
1-4	MODIF.KIT INCU.8000 SC, USA	2M20664	
1-2	DEPOSIT TRAY,COMPL.		2M20320
1	DEPOSIT PLATE	2M20695	
2	RUBBER BUFFER	2M20704	
3	LABEL	2M20631	
4	SET OF LABELS USA		2M20580

## MODIFICATION KIT CORNER PIECE

Picture 72



Item No.	Description	Part No.	Order-Code
1-10	MODIFICATION KIT CORNER PIECE		2M20722
1	CORNER, LEFT	2M20650	
2	ANGLE	2M20683	
3	SPRING	D14967	
4	GROOVED PIN DIN 1473-8X25-A2	1342053	
5	CORNER, RIGHT	2M20651	
6	MOUNTING INSTRUCTIONS, without illustration	2M20721	
8	INSERT, without illustration	8711067	
9	STIFFENING		2M20214
10	SCREW DWN 562-3X12-F-A2-H	1338986	